
Getting upstairs in safety

Stair Lift



RL-50 – The space – saving Wheelchair Ceiling Lift System

The Högg Wheelchair Ceiling Lift is a unique solution for the wheelchair user to be completely independent and mobile around the house. It allows the user to remain within their existing home rather than face the stress of moving house. The sturdy quality Swiss made construction, is in accordance with the highest safety standards to ensure comfortable and safe transportation between floors, indeed you are not limited to just going up only one floor. Boarding, departure and parking points can be adapted to suit your needs and local conditions. Attaching the wheelchair to the ceiling lift is simple and quick, requiring little physical effort, you are then ready to go.









Parking position of lift.

Safe and simple attachment of wheelchair.

Comfortable driving position.

The Högg RL-50 can be installed almost anywhere, even in narrow steep and curved stairways, with some adjustments it can even be made to go through doorways. The drive mechanism employs the latest technology for your guarantee of long-term safety.

Because the RL-50 uses the ceiling lift system developed by Högg of Switzerland, the stairs and passages remain free of obstruction, creating a safe and unobtrusive layout.

The RL-50 can also be equipped with a transportation carriage as an optional extra. This can be used for the carriage of shopping bags and other household items.







Technical details

Manufacturer: Högg Liftsysteme AG, Switzerland

Nominal load: 120 kg

Controls: key switch, Hold-on scanner (pressure scanner)

Emergency stop

Drive system: 42 V AC/pinion gear on chain

Rail: steel rail 35×40 mm

Electrical connection: 230 V AC

Fuse 13 A

Speed: 0,1 m/sec

Technical conformance: CE/Technical Inspection Authority

approved

Area of application: inclines from 0 degrees (horizontal)

to 60 degrees

Standard colour: silve

Standard design: includes wheelchair


attachment bracket


Optional: special mounting

Special seat with wheels Transport carriage

Cordless call and transmission points

Subject to changes without notice.





All measurements in mm.



Bürgistrasse 15 CH-9620 Lichtensteig Phone 071-987 66 80 Fax 071-987 66 89 lift@hoegg.ch www.hoegg.ch